

# Anti-DHRS9 antibody



Catalog Number: 142235

## Product name

Anti-DHRS9 antibody

## Immunogen

[Human DHRS9/3-alpha-HSD \(His Tag\) recombinant protein](#)

## Specificity

Human DHRS9 / 3-alpha-HSD

## Antibody description

Rabbit Polyclonal to Human DHRS9

## Preparation

Produced in rabbits immunized with purified, recombinant Human DHRS9 / 3-alpha-HSD . DHRS9 / 3-alpha-HSD specific IgG was purified by Human DHRS9 / 3-alpha-HSD affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, ELISA, IP

## Dilutions


WB: 10-20 µg/ml

ELISA: 0.1-0.2 µg/ml


This antibody can be used at 0.1-0.2 µg/ml with the appropriate secondary reagents to detect Human DHRS9 / 3-alpha-HSD.

IP: 1-4 uL/mg of lysate

## Validations



DHRS9 / 3-alpha-HSD Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot



DHRS9 / 3-alpha-HSD Antibody, Rabbit PAb, Antigen Affinity Purified, Immunoprecipitation