

Anti-MMP9 antibody



Catalog Number: 142247

Product name

Anti-MMP9 antibody

Immunogen

[Rat MMP-9/MMP9 \(His Tag\) recombinant protein](#)

Specificity

Rat MMP-9 / MMP9

Antibody description

Rabbit Polyclonal to Rat MMP9

Preparation

Produced in rabbits immunized with purified, recombinant Rat MMP-9 / MMP9 . MMP-9 / MMP9 specific IgG was purified by Rat MMP-9 / MMP9 affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC-P


Dilutions

ELISA: 0.1-0.2 µg/ml

This antibody can be used at 0.1-0.2 µg/ml with the appropriate secondary reagents to detect Rat MMP-9 / MMP9.

IHC-P: 0.1-2 µg/mL

Validations



MMP-9 / MMP9 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry