

Anti-HA antibody



Catalog Number: 142286

Product name

Anti-HA antibody

Immunogen

[H5N1 HA \(His Tag\) recombinant protein](#)

Specificity

Influenza A H5N1 (A/turkey/Turkey/1/2005)
Hemagglutinin / HA

Antibody description

Rabbit Polyclonal to H5N1 HA

Preparation

Produced in rabbits immunized with purified, recombinant Influenza A H5N1 (A/turkey/Turkey/1/2005) Hemagglutinin / HA . Influenza A H5N1 (A/turkey/Turkey/1/2005) Hemagglutinin / HA specific IgG was purified by Influenza A H5N1 (A/turkey/Turkey/1/2005) Hemagglutinin / HA affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ELISA, ICC/IF, IF, IP


Dilutions

WB: 0.5-5 µg/ml

ELISA: 0.1-0.2 µg/ml

This antibody can be used at 0.1-0.2 µg/ml with the appropriate secondary reagents to detect Influenza A H5N1 (A/turkey/Turkey/1/2005) Hemagglutinin / HA.

Validations



Lanes	A	B
Items		
Sample (Recombinant protein)	Turkey-HA-His (Cat#11061-V08H1)	
Sample Volume (ng/lane)	50	10
Gel	13% SDS-PAGE reducing gel	
Recommended Concentration	0.5-5 µg/ml	
Secondary Antibody	Dylight 800-labeled Antibody to Rabbit IgG (1+1), at 1:5000 dilution.	
Developed using Odyssey imaging system.		

Influenza A H5N1 (A/turkey/Turkey/1/2005) Hemagglutinin / HA Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot