

Anti-HA antibody

Catalog Number: 142306

Product name

Anti-HA antibody

Immunogen

[H5N1 HA \(His Tag\) recombinant protein](#)

Specificity

Influenza A H5N1 (A/Vietnam/1194/2004)
Hemagglutinin / HA

Antibody description

Rabbit Polyclonal to H5N1 HA

Preparation

Produced in rabbits immunized with purified, recombinant Influenza A H5N1 (A/Vietnam/1194/2004) Hemagglutinin / HA . Influenza A H5N1 (A/Vietnam/1194/2004) Hemagglutinin / HA specific IgG was purified by Influenza A H5N1 (A/Vietnam/1194/2004) Hemagglutinin / HA affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ELISA, ICC/IF, IF, IP


Dilutions

WB: 0.5-5 µg/ml

ELISA: 0.1-0.2 µg/ml

This antibody can be used at 0.1-0.2 µg/ml with the appropriate secondary reagents to detect Influenza A H5N1 (A/Vietnam/1194/2004) Hemagglutinin / HA.

Validations



Items	Lanes	A	B
Sample (Recombinant protein)	Vietnam-HA-His (Cat#11062-V08H1)		
Sample Volume (ng/lane)	50	10	
Gel	13% SDS-PAGE reducing gel		
Recommended Concentration	0.5-5 µg/ml		
Secondary Antibody	Dylight 800-labeled Antibody To Rabbit IgG (H+L), at 1:5000 dilution.		
Developed using Odyssey imaging system.			

Influenza A H5N1 (A/Vietnam/1194/2004) Hemagglutinin / HA Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot