

# Anti-STK16 antibody



Catalog Number: 142310

## Product name

Anti-STK16 antibody

## Immunogen

[Human STK16/PKL12 \(His & NusA Tag\) recombinant protein](#)

## Specificity

Human STK16 / PKL12

## Antibody description

Rabbit Polyclonal to Human STK16

## Preparation

Produced in rabbits immunized with purified, recombinant Human STK16 / PKL12 . STK16 / PKL12 specific IgG was purified by Human STK16 / PKL12 affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IHC-P

## Dilutions

IHC-P: 0.1-2 µg/mL

## Validations



Human STK16/PKL12 Immunohistochemistry(IHC) 14712