

Anti-SAA4 antibody



Catalog Number: 142327

Product name

Anti-SAA4 antibody

Immunogen

[Human SAA4 \(GST Tag\) recombinant protein](#)

Specificity

Human SAA4

Antibody description

Mouse Monoclonal to Human SAA4

Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human SAA4 and conjugated with PE under optimum conditions, the unreacted PE was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from

prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

FCM

Dilutions

FCM

Validations



SAA4 Antibody, Mouse MAb, Flow Cytometry