

# Anti-CD33 antibody



Catalog Number: 142454

## Product name

Anti-CD33 antibody

## Immunogen

[Human CD33 / Siglec-3 \(His Tag\) recombinant protein](#)

## Specificity

Human CD33 / Siglec-3

## Antibody description

Mouse Monoclonal to Human CD33

## Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human CD33 / Siglec-3 and conjugated with PerCP under optimum conditions, the unreacted PerCP was removed.

## Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

## Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from

prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1


## Applications

FCM

## Dilutions

FCM

## Validations



CD33 / Siglec-3 Antibody (PerCP), Mouse MAb, Flow Cytometry