

# Anti-PECAM1 antibody



Catalog Number: 142466

## Product name

Anti-PECAM1 antibody

## Immunogen

[Human CD31/PECAM-1 \(His Tag\) recombinant protein](#)

## Specificity

Human CD31 / PECAM-1

## Antibody description

Rabbit Monoclonal to Human PECAM1

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human CD31 / PECAM-1 .

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid

contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

ICC/IF, IF

## Dilutions

ICC/IF: 10-25 µg/mL

## Validations



CD31 / PECAM-1 Antibody, Rabbit MAb, Immunofluorescence