

Anti-HA antibody



Catalog Number: 142507

Product name

Anti-HA antibody

Immunogen

[Influenza B HA \(His Tag\) recombinant protein](#)

Specificity

Influenza B (B/Brisbane/60/2008) Hemagglutinin / HA1

Antibody description

Rabbit Polyclonal to Influenza B HA

Preparation

Produced in rabbits immunized with purified, recombinant Influenza B (B/Brisbane/60/2008) Hemagglutinin / HA1 . Influenza B (B/Brisbane/60/2008) Hemagglutinin / HA1 specific IgG was purified by Influenza B (B/Brisbane/60/2008) Hemagglutinin / HA1 affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB, ICC/IF, IF, IP

Dilutions

WB: 0.5-5 µg/ml

Validations



Items	Lanes	A	B
Sample (Recombinant protein)		B-B5-08-HA1-His (Cat#40016-V08H1)	
Sample Volume (ng/lane)		50	10
Gel		13% SDS-PAGE reducing gel	
Recommended Concentration		0.5-5 µg/ml	
Secondary Antibody		Dylight 800-labeled Antibody To Rabbit IgG (H+L), at 1:5000 dilution. Developed using Odyssey imaging system.	

Influenza B (B/Brisbane/60/2008) Hemagglutinin / HA1 Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot