

# Anti-HA antibody



Catalog Number: 142560

## Product name

Anti-HA antibody

## Immunogen

[H3N2 HA \(His Tag\) recombinant protein](#)

## Specificity

Influenza A H3N2 (A/Aichi/2/1968) Hemagglutinin / HA1

## Antibody description

Rabbit Polyclonal to H3N2 HA

## Preparation

Produced in rabbits immunized with purified, recombinant Influenza A H3N2 (A/Aichi/2/1968) Hemagglutinin / HA1 . Influenza A H3N2 (A/Aichi/2/1968) Hemagglutinin / HA1 specific IgG was purified by Influenza A H3N2 (A/Aichi/2/1968) Hemagglutinin / HA1 affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB, ICC/IF, IF, IP

## Dilutions

WB: 0.5-5 µg/ml

## Validations



Items	Lanes	A	B
Sample (Recombinant protein)		H3N2-AI-68-HA1-His (Cat#11707-V08H1)	
Sample Volume (ng/lane)		50	10
Gel		1.3% SDS-PAGE reducing gel	
Recommended Concentration		0.5-5 µg/ml	
Secondary Antibody		Dylight 800-labeled Antibody To Rabbit IgG (H+L), at 1:5000 dilution.	
Developed using Odyssey imaging system.			

Influenza A H3N2 (A/Aichi/2/1968) Hemagglutinin / HA1 Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot