

# Anti-FBXO11 antibody



Catalog Number: 142563

## Product name

Anti-FBXO11 antibody

## Specificity

Human FBXO11

## Antibody description

Rabbit Polyclonal to Human FBXO11

## Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the N-terminus of the Human FBXO11, and purified by antigen affinity chromatography.

## Formulation

0.2  $\mu$ m filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It

is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 5-20  $\mu$ g/ml

## Validations



Lanes	A	B
Sample (whole cell lysate)	293T	HeLa
Sample Volume ( $\mu$ g/lane)	30	30
Gel	7.5%SDS-PAGE reducing gel	
Recommended Concentration	5-20 $\mu$ g/ml	
Secondary Antibody	Dylight 800 labeled Antibody to Rabbit IgG (H+), at 1:5000 dilution.	
Explanation	Developed using Odyssey imaging system. Predicted band size : 94 kDa Observed band size : 120 kDa Additional bands at : 58 kDa, 80 kDa and 45 kDa (We are unsure as to the identity of these extra bands.)	

Human FBXO11 Western blot (WB) 14628