

Anti-PDCD4 antibody



Catalog Number: 142581

Product name

Anti-PDCD4 antibody

Specificity

Human PDCD4

Antibody description

Rabbit Polyclonal to Human PDCD4

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the center region of the Human PDCD4, and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid

contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ICC/IF, IF

Dilutions

ICC/IF: 0.5-1.5 µg/mL

Validations



PDCD4 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunofluorescence