

Anti-HNF4A antibody



Catalog Number: 142600

Product name

Anti-HNF4A antibody

Specificity

Human HNF4A

Antibody description

Rabbit Polyclonal to Human HNF4A

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the N-terminus of the Human HNF4A, and purified by antigen affinity chromatography.

Formulation

0.2 μ m filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, ICC/IF, IF

Dilutions

WB: 5-20 μ g/ml


ICC/IF: 0.5-3 μ g/mL

Validations



Lanes	A
Items	
Sample (whole cell lysate)	A549
Sample Volume (μ g/lane)	30
Gel	13% SDS-PAGE reducing gel
Recommended Concentration	5-20 μ g/ml
Secondary Antibody	Dylight 800-labeled Antibody to Rabbit IgG (H+L), at 1:5000 dilution.
Explanation	Developed using Odyssey imaging system. Predicted band size: 53 kDa Observed band size: 57 kDa Additional bands at: 65 kDa We are unsure as to the identity of these extra bands.

HNF4A Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot



HNF4A Antibody, Rabbit PAb, Antigen Affinity Purified, Immunofluorescence