

Anti-HA antibody



Catalog Number: 142626

Product name

Anti-HA antibody

Immunogen

[Influenza B HA \(His Tag\) recombinant protein](#)

Specificity

Influenza B (B/Massachusetts/03/2010)
Hemagglutinin / HA1

Antibody description

Rabbit Polyclonal to Influenza B HA

Preparation

Produced in rabbits immunized with purified, recombinant Influenza B (B/Massachusetts/03/2010) Hemagglutinin / HA1 . Influenza B (B/Massachusetts/03/2010) Hemagglutinin / HA1 specific IgG was purified by Influenza B (B/Massachusetts/03/2010) Hemagglutinin / HA1 affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ELISA, IHC-P


Dilutions

WB: 0.5-1 µg/ml

ELISA: 0.1-0.2 µg/ml

This antibody can be used at 0.1-0.2 µg/ml with the appropriate secondary reagents to detect Influenza B (B/Massachusetts/03/2010) Hemagglutinin / HA1.

Validations



Items	Lanes	A	B	C
Sample (Recombinant protein)	B-MC-10 HA1-His (Cat#40191-V08H1)			
Sample Volume (µg/lane)		20	5	1
Gel	13% SDS-PAGE reducing gel			
Recommended Concentration	0.5-1 µg/ml			
Secondary Antibody	Dylight 800-labeled Antibody to Rabbit IgG (H+L), at 1:5000 dilution.			
Developed using Odyssey imaging system.				

Influenza B (B/Massachusetts/03/2010) Hemagglutinin / HA1 Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot