

Anti-MAF antibody



Catalog Number: 142697

Product name

Anti-MAF antibody

Specificity

Human MAF

Antibody description

Rabbit Polyclonal to Human MAF

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the N-terminus of the Human MAF, and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It

is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IP

Dilutions

IP: 1-4 µl/mg of lysate

Validations



Item	Lane	A
Sample (wholecell lysate)		293T
Sample quantity		0.5 mg
IP antibody quantity		4 µl
Protein G agarose		15 µl of 50% Protein G agarose
Gel		13% SDS-PAGE reducing gel
Primary antibody		KLJ1-SMCC-SINOAA2568 KLJ1-SMCC-SINOAA2788 antibody at 10 µg/ml
Secondary antibody		DyLight 800 labeled antibody to rabbit IgG (H+L), at 1:5000 dilution

Human MAF Immunoprecipitation(IP) 14673