

Anti-HA antibody



Catalog Number: 142702

Product name

Anti-HA antibody

Immunogen

[H5N1 HA \(28 Ser/Trp, His Tag\) recombinant protein](#)

Specificity

Influenza H5N1 (A/Hong kong/213/2003)
Hemagglutinin / HA

Antibody description

Rabbit Polyclonal to H5N1 HA

Preparation

Produced in rabbits immunized with purified, recombinant Influenza H5N1 (A/Hong kong/213/2003) Hemagglutinin / HA . Influenza H5N1 (A/Hong kong/213/2003) Hemagglutinin / HA specific IgG was purified by Influenza H5N1 (A/Hong kong/213/2003) Hemagglutinin / HA affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It

is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ELISA, FCM, ICC/IF, IF, IP


Dilutions

WB: 0.5-5 µg/ml

ELISA: 0.1-0.2 µg/ml

This antibody can be used at 0.1-0.2 µg/ml with the appropriate secondary reagents to detect Influenza H5N1 (A/Hong kong/213/2003) Hemagglutinin / HA.

Validations



Items	Lanes	A	B
Sample (Recombinant protein)		H5N1-HK-03-HA-His(E) (Cat#11713-V08H)	
Sample Volume (ng/lane)		50	10
Gel		13% SDS-PAGE reducing gel	
Recommended Concentration		0.5-5 µg/ml	
Secondary Antibody		Dylight 800-labeled Antibody to Rabbit IgG (H+L), at 1:5000 dilution. <small>Developed using Odyssey Imaging system.</small>	

Influenza H5N1 (A/Hong kong/213/2003) Hemagglutinin / HA Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot