

# Anti-ACSS3 antibody



Catalog Number: 142708

## Product name

Anti-ACSS3 antibody

## Specificity

Human ACSS3

## Antibody description

Rabbit Polyclonal to Human ACSS3

## Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the center region of the Human ACSS3, and purified by antigen affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


IHC-P, IP

## Dilutions


IHC-P: 0.1-2 µg/mL

IP: 1-4 µl/mg of lysate

## Validations



## Human ACSS3 Immunohistochemistry(IHC) 14677



Items	Lanes	A
Sample (whole cell lysate)		HepG2
Sample quantity		0.5 mg
IP antibody quantity		4 µL
Protein G agarose		15 µl of 50% Protein G Agarose
Gel		7.5% SDS-PAGE/reducing gel
Primary antibody		KLH-SMCC-SINO2632/KLH-SMCC-SINO2639 antibody at 10 µg/ml
Secondary antibody		DyLight 800-labeled antibody to rabbit IgG (H+L) at 1:5000 dilution

## Human ACSS3 Immunoprecipitation(IP) 14678