

Anti-IL1A antibody



Catalog Number: 142854

Product name

Anti-IL1A antibody

Immunogen

[Mouse IL-1 alpha / IL-1F1 / IL1A recombinant protein](#)

Specificity

Mouse IL1A/IL-1A/IL-1F1

Antibody description

Rabbit Polyclonal to Mouse IL1A

Preparation

Produced in rabbits immunized with purified, recombinant Mouse IL1A/IL-1A/IL-1F1 . IL1A/IL-1A/IL-1F1 specific IgG was purified by Mouse IL1A/IL-1A/IL-1F1 affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ELISA, IP

Dilutions


WB: 5-20 µg/ml

ELISA: 0.1-0.2 µg/mL


This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Mouse IL1A/IL-1A/IL-1F1.

IP: 4-8 uL/mg of lysate

Validations



Mouse IL1A/IL-1A/IL-1F1 Western blot (WB)



15852

Mouse IL1A/IL-1A/IL-1F1 Immunoprecipitation(IP)
15853