

iNOS Polyclonal Antibody



Catalog Number: 160058

Product name

Rabbit IgG

iNOS Polyclonal Antibody

Applications

Specificity

WB

Human, Mouse, Rat

Dilutions

Antibody description


WB 1:500 - 1:2000

Polyclonal antibody to iNOS

Validations

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 500-730 of human iNOS (NP_000616.3).



Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Western blot - iNOS Polyclonal Antibody

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Western blot analysis of extracts of Raw 264.7 cells, using iNOS antibody at 1:1000 dilution. Raw264.7 cells were treated by LPS (1 µg/ml) for 8 hours. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.

Clonality

Polyclonal

Ig Type