

# TLR2 Polyclonal Antibody

Catalog Number: 160065

## Product name

TLR2 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to TLR2

## Preparation

Antigen: Recombinant protein of human TLR2

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

## Validations




Western blot - TLR2 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using TLR2 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunohistochemistry - TLR2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat lung using TLR2 antibody at dilution of 1:200 (40x



lens).

Immunohistochemistry - TLR2 Polyclonal Antibody

# TLR2 Polyclonal Antibody



Catalog Number: 160065

---

Immunohistochemistry of paraffin-embedded rat

bone marrow using TLR2 antibody at dilution of  
1:200 (40x lens).