

ATXN1 Polyclonal Antibody



Catalog Number: 160095

Product name

ATXN1 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to ATXN1

Preparation

Antigen: A synthetic peptide of human ATXN1

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:200 - 1:500

IF 1:20 - 1:50

Validations



Western blot - ATXN1 Polyclonal Antibody

Western blot analysis of extracts of CEM cells, using ATXN1 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.