

# RIG-I / DDX58 Polyclonal Antibody



Catalog Number: 160101

## Product name

RIG-I / DDX58 Polyclonal Antibody

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to RIG-I / DDX58

## Preparation

Antigen: A synthetic peptide of human RIG-I / DDX58

## Formulation

PBS with 0.02% sodium azide, pH7.3.

## Storage

Store at 4°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:1000

## Validations



Western blot - RIG-I / DDX58 Polyclonal Antibody

Western blot analysis of extracts of HUVEC cells, using RIG-I / DDX58 antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.