

ELL Polyclonal Antibody



Catalog Number: 160126

Product name

ELL Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to ELL

Preparation

Antigen: A synthetic peptide of human ELL

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:200 - 1:500

Validations



Western blot - ELL Polyclonal Antibody

Western blot analysis of extracts of 293 cells, using ELL antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.