

CASP7 Polyclonal Antibody



Catalog Number: 160140

Product name

CASP7 Polyclonal Antibody

Specificity

Human, Rat

Antibody description

Polyclonal antibody to CASP7

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 24-303 of human CASP7 (NP_001218.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


IHC, IP

Dilutions

IHC 1:50 - 1:100

IP 1:50 - 1:100

Validations



Immunohistochemistry - CASP7 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat heart using CASP7 antibody at dilution of 1:100 (40x lens).