

APAF1 Polyclonal Antibody



Catalog Number: 160148

Product name

APAF1 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to APAF1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human APAF1 (NP_863651.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - APAF1 Polyclonal Antibody

Western blot analysis of extracts of 293 cells, using APAF1 antibody .