

# ADAM17 Polyclonal Antibody



Catalog Number: 160159

## Product name

ADAM17 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to ADAM17

## Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 700 to the C-terminus of human ADAM17 (NP\_003174.3).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

WB 1:500 - 1:2000

IHC 1:100 - 1:200


IF 1:50 - 1:200

## Validations



Western blot - ADAM17 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using ADAM17 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 5s.



Immunohistochemistry - ADAM17 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat heart using ADAM17 antibody at dilution of 1:100 (40x lens).