

BIRC3 Polyclonal Antibody



Catalog Number: 160161

Product name

BIRC3 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to BIRC3

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human BIRC3 (NP_001156.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:200 - 1:2000


IF 1:50 - 1:200

Validations



Western blot - BIRC3 Polyclonal Antibody

Western blot analysis of extracts of U-937 cells, using BIRC3 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunofluorescence - BIRC3 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using BIRC3 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.