

BIRC2 Polyclonal Antibody



Catalog Number: 160168

Product name

BIRC2 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to BIRC2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human BIRC2 (NP_001157.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:200

Validations



Western blot - BIRC2 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using BIRC2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.



Immunofluorescence - BIRC2 Polyclonal Antibody

Immunofluorescence analysis of A549 cells using BIRC2 antibody.