

LRRK2 Polyclonal Antibody



Catalog Number: 160169

Product name

LRRK2 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to LRRK2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 2355-2527 of human LRRK2 (NP_940980.3).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:200

Validations



Immunofluorescence - LRRK2 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using LRRK2 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.