

TRAF3 Polyclonal Antibody



Catalog Number: 160171

Product name

TRAF3 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to TRAF3

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human TRAF3 (NP_003291.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:200 - 1:1000


IHC 1:50 - 1:200

Validations



Western blot - TRAF3 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using TRAF3 antibody at 1:500 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunohistochemistry - TRAF3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using TRAF3 antibody at dilution of 1:100 (40x lens).