

# CAPN1 Polyclonal Antibody



Catalog Number: 160181

## Product name

CAPN1 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to CAPN1

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 20-330 of human CAPN1 (NP\_005177.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:2000

## Validations



Western blot - CAPN1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using CAPN1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.