

CCNT1 Polyclonal Antibody



Catalog Number: 160190

Product name

CCNT1 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to CCNT1

Preparation

Antigen: A synthetic peptide of human CCNT1

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, FC

Dilutions

WB 1:200 - 1:500

FC 1:20 - 1:50

Validations



Western blot - CCNT1 Polyclonal Antibody

Western blot analysis of extracts of K-562 cells, using CCNT1 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.