

# NCK1 Polyclonal Antibody



Catalog Number: 160207

## Product name

NCK1 Polyclonal Antibody

WB 1:500 - 1:2000

IHC 1:50 - 1:100

IF 1:50 - 1:100

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to NCK1

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 58-377 of human NCK1 (NP\_006144.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type


Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

## Validations



Western blot - NCK1 Polyclonal Antibody

Western blot analysis of extracts of mouse brain, using NCK1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time:



Immunofluorescence - NCK1 Polyclonal Antibody

Immunofluorescence analysis of A549 cells using NCK1 antibody. Blue: DAPI for nuclear staining.