

MUC16 Polyclonal Antibody



Catalog Number: 160211

Product name

MUC16 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to MUC16

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 13965-14064 of human MUC16 (NP_078966.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC

Dilutions


IHC 1:50 - 1:200

Validations



Immunohistochemistry - MUC16 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse lung using MUC16 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - MUC16 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver injury using MUC16 antibody at dilution of 1:100 (40x lens).