

# NPC1L1 Polyclonal Antibody



Catalog Number: 160240

## Product name

NPC1L1 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to NPC1L1

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1030-1100 of human NPC1L1 (NP\_001095118.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:1000 - 1:4000

## Validations



## Western blot - NPC1L1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using NPC1L1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 15s.