

FCGR2C Polyclonal Antibody



Catalog Number: 160254

Product name

Rabbit IgG

FCGR2C Polyclonal Antibody

Applications

Specificity

WB

Human

Dilutions

Antibody description

WB 1:1000 - 1:4000

Polyclonal antibody to FCGR2C

Validations

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 43-223 of human FCGR2C (NP_963857.3).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type



Western blot - FCGR2C Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using FCGR2C antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.