

ATP2B4 Polyclonal Antibody



Catalog Number: 160282

Product name

ATP2B4 Polyclonal Antibody

Specificity

Mouse

Antibody description

Polyclonal antibody to ATP2B4

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1141-1205 of human ATP2B4 (NP_001675.3).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB

Dilutions

WB 1:1000 - 1:4000

Validations

Mouse brain



Western blot - ATP2B4 Polyclonal Antibody

Western blot analysis of extracts of mouse brain, using ATP2B4 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5s.