

CECR2 Polyclonal Antibody



Catalog Number: 160314

Product name

Rabbit IgG

CECR2 Polyclonal Antibody

Applications

Specificity

WB

Mouse, Rat

Dilutions

Antibody description

WB 1:200 - 1:1000

Polyclonal antibody to CECR2

Validations

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 170-370 of human CECR2 (NP_113601.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type



Western blot - CECR2 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using CECR2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 60s.