

# PPAR alpha Monoclonal Antibody



Catalog Number: 160336

## Product name

PPAR alpha Monoclonal Antibody

## Specificity

Human, Mouse

## Antibody description

Monoclonal antibody to PPAR alpha

## Preparation

Antigen: Recombinant protein of human PPAR alpha

## Formulation

PBS with 0.02% sodium azide, pH7.3.

## Storage

Store at 4°C. Avoid freeze / thaw cycles.

## Clonality

Monoclonal

## Ig Type

Mouse IgG

## Applications

WB, IHC, IF, FC

## Dilutions


WB 1:500 - 1:1000

IHC 1:10 - 1:25

IF 1:10 - 1:25

FC 1:10 - 1:25

## Validations



Western blot - PPAR alpha Monoclonal Antibody

Western blot analysis of extracts of various cell lines, using PPAR alpha antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.