

MS4A8B Polyclonal Antibody



Catalog Number: 160424

Product name

MS4A8B Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to MS4A8B

Preparation

Antigen: A synthetic peptide of human MS4A8B

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations



Western blot - MS4A8B Polyclonal Antibody

Western blot analysis of extracts of U-87MG cells, using MS4A8B antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.



Immunohistochemistry - MS4A8B Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse heart using MS4A8B antibody at dilution of 1:100 (40x lens).