

TMEM74 Polyclonal Antibody



Catalog Number: 160458

Product name

TMEM74 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to TMEM74

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human TMEM74 (NP_694560.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - TMEM74 Polyclonal Antibody

Western blot analysis of extracts of DU145 cells, using TMEM74 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.