

# VDAC3 Polyclonal Antibody



Catalog Number: 160546

## Product name

VDAC3 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to VDAC3

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-283 of human VDAC3 (NP\_005653.3).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:1000 - 1:2000


IF 1:50 - 1:200

## Validations



### Western blot - VDAC3 Polyclonal Antibody

Western blot analysis of extracts of A-431 cells, using VDAC3 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5s.



### Immunofluorescence - VDAC3 Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using VDAC3 antibody at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.