

LYZL6 Polyclonal Antibody



Catalog Number: 160588

Product name

LYZL6 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to LYZL6

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 20-148 of human LYZL6 (NP_065159.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:200

Validations



Western blot - LYZL6 Polyclonal Antibody

Western blot analysis of extracts of human testis, using LYZL6 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunofluorescence - LYZL6 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using LYZL6 antibody. Blue: DAPI for nuclear staining.