

CYBA Polyclonal Antibody



Catalog Number: 160590

Product name

CYBA Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to CYBA

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 126-195 of human CYBA (NP_000092.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:50 - 1:100

Validations



Immunofluorescence - CYBA Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using CYBA antibody at dilution of 1:100 (40x lens).
Blue: DAPI for nuclear staining.