

# HNRNPDL Polyclonal Antibody



Catalog Number: 160594

## Product name

HNRNPDL Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to HNRNPDL

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 261-420 of human HNRNPDL (NP\_112740.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:100

IF 1:50 - 1:100

## Validations



Western blot - HNRNPDL Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using HNRNPDL antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL

Basic Kit .Exposure time: 30s.

Immunofluorescence - HNRNPDL Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using HNRNPDL antibody at dilution of 1:100. Blue: DAPI for nuclear staining.

