

DEFB132 Polyclonal Antibody



Catalog Number: 160597

Product name

DEFB132 Polyclonal Antibody

Specificity

Mouse

Antibody description

Polyclonal antibody to DEFB132

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-95 of human DEFB132 (NP_997352.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - DEFB132 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using DEFB132 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.