

KLK14 Polyclonal Antibody



Catalog Number: 160608

Product name

KLK14 Polyclonal Antibody

Specificity

Mouse

Antibody description

Polyclonal antibody to KLK14

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 41-267 of human KLK14 (NP_071329.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB

Dilutions

WB 1:500 - 1:2000

Validations

Mouse stomach



Western blot - KLK14 Polyclonal Antibody

Western blot analysis of extracts of mouse stomach, using KLK14 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.