

# NUAK2 Polyclonal Antibody



Catalog Number: 160618

## Product name

NUAK2 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to NUAK2

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 523-672 of human NUAK2 (NP\_112214.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:2000

## Validations



Western blot - NUAK2 Polyclonal Antibody

Western blot analysis of extracts of rat liver, using NUAK2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.