

COL17A1 Monoclonal Antibody



Catalog Number: 160630

Product name

COL17A1 Monoclonal Antibody

IHC 1:50 - 1:200

IF 1:50 - 1:200

IP 1:20 - 1:50

Specificity

Human, Mouse, Rat

Antibody description

Monoclonal antibody to COL17A1

Preparation

Antigen: Recombinant protein of human COL17A1

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

WB, IHC, IF, IP

Dilutions


WB 1:500 - 1:2000

Validations



Immunohistochemistry - COL17A1 Monoclonal Antibody

Immunohistochemistry of paraffin-embedded mouse skin using COL17A1 antibody at dilution of 1:100 (40x lens).



Immunofluorescence - COL17A1 Monoclonal Antibody

Immunofluorescence analysis of HUVEC cells using COL17A1 antibody .