

BCL2L13 Polyclonal Antibody



Catalog Number: 160687

Product name

BCL2L13 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to BCL2L13

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human BCL2L13 (NP_056182.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000


IF 1:20 - 1:50

Validations



Western blot - BCL2L13 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using BCL2L13 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunofluorescence - BCL2L13 Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using BCL2L13 antibody. Blue: DAPI for nuclear staining.